

S-Path **NEW** BILIARY STENT PLACEMENT SOLUTION

With Novel Integrated **RX-OTW** Self-Detaching Delivery System

Pre-mounted single stent insertion via **one system pass**

- Secures accurate guidewire access
without any manipulation
- *Reduces procedure steps and time!*



Endo GI Medical
ADVANCED GASTROINTESTINAL SOLUTIONS

Drawbacks of Current Biliary Stent Delivery Systems


Lost guidewires

Multiple entry points for multiple stents

Prolonged Procedure time

Procedural complications

Discover The S-Path to Improved Procedural Success with The EndoGI Superior Technology Platform



S-PATH STENT TECHNOLOGY FEATURE

Novel integrated RX OTW self-detaching delivery system

No guidewire maneuvering
Uniquely designed tip

Guidewire compatibility

Endoscopically visible radiopaque polymer

CLINICAL & PROCEDURAL BENEFIT

Secures guidewire access
Accurate guidewire positioning

Reduced tissue trauma with better pushability,
ductal access and surgical outcome

Compatible with all 0.035“(0.89 mm) guidewires
and 10 FR 8.5 FR & 7 FR diameters
Consistent with existing surgical procedures &
protocols

Good procedural visibility



www.endo-gi.com

EndoGI Medical 12, Wadi El Headj Nazareth Zip 1603611